

FORSURE

BLOOD GROUPING REAGENTS


FORSURE Blood Grouping Reagents are used for the identification of blood types. The test procedure is based on agglutination principle. Our blood grouping reagents are prepared from in-vitro culture supernatants of hybridized immunoglobulin secreting mouse cell lines. Anti-A, Anti-B and Anti-AB reagents (here under referred to as ABO reagents) are used for qualitative in-vitro determination of human blood groups of the ABO system to determine the blood type. Anti-D reagent is used for the qualitative determination of Rhesus factor on human blood groups. These reagents can be used in slide, tube and microplate methods.

Our blood grouping reagents range also include the Enhancement Media such as Anti-Human Globulin, Bovine Albumin, Negative Control for Anti-D Reagent and Low Ionic Strength Saline (LISS).

FORSURE

BLOOD GROUPING REAGENTS

Product Code	Blood Grouping Reagents	Size
FS-A256-B10-010	Forsure Anti-A Monoclonal Reagent (1/256), 10m/vial	10 vials/pack
FS-A512-B10-010	Forsure Anti-A Monoclonal Reagent (1/512), 10ml/vial	10 vials/pack
FS-B256-B10-010	Forsure Anti-B Monoclonal Reagent (1/256), 10ml/vial	10 vials/pack
FS-B512-B10-010	Forsure Anti-B Monoclonal Reagent (1/512), 10ml/vial	10 vials/pack
FS-AB256-B10-010	Forsure Anti-AB Monoclonal Reagent (1/256), 10ml/vial	10 vials/pack
FS-AB512-B10-010	Forsure Anti-AB Monoclonal Reagent (1/512), 10ml/vial	10 vials/pack
FS-D64-B10-010	Forsure Anti-D IgG/IgM Blend Reagent (1/64), 10ml/vial	10 vials/pack
FS-D128-B10-010	Forsure Anti-D IgG/IgM Blend Reagent (1/128), 10ml/vial	10 vials/pack
FS-BA22-B10-001	Forsure Bovine Albumin 22%, 10ml/vial	1 vial/box
FS-BA30-B10-001	Forsure Bovine Albumin 30%, 10ml/vial	1 vial/box
FS-G256-B10-001	Forsure Anti-Human Globulin (1/256), 10ml/vial	1 vial/box
FS-RhC-B10-001	Forsure Negative Control for Anti-D Reagent, 10ml/vial	1 vial/box
FS-LISS-B10-001	Forsure Low Ionic Strength Saline, 10ml/vial	1 vial/box
FS-ABO-B10-003	Forsure ABO Set [Anti-A (1/256), Anti-B (1/256), Anti-D (1/64)], 10ml/vial	3 vials/box
FS-ABO-B10-004	Forsure ABO Set [Anti-A (1/512), Anti-B (1/512), Anti-AB (1/512), Anti-D (1/128)], 10ml/vial	4 vials/box



Forsure ABO Set



** single vial packing

TS Diagnostics Sdn Bhd

No 7 & 9, Jalan TPP 5/11,
Taman Perindustrian Puchong,
47160 Puchong, Selangor,
Malaysia.

Ph number +6 03 8068 7311
Fax number +6 03 8051 1778
contactus@forsure-test.com
www.forsure-test.com

