



FORSURE

FERTILITY RAPID TEST

FORSURE Fertility Rapid Test range are based on lateral flow immunoassay for the detection of hormones including Human chorionic gonadotropin (hCG) for pregnancy, Luteinizing hormone (LH) for ovulation and Follicle-stimulating hormone (FSH) for menopause in urine samples. These tests come in various formats and are conveniently packed in aluminium pouch to suit the lab, point of care as well as home use.

FORSURE Fertility Rapid Test range are accurate, reliable, easy to use, competitively priced and results are obtained between 1 to 5 minutes.

forsure

FERTILITY RAPID TEST

Product Code	Description	Format	Size
FS-HCG-S30-001b	Forsure Pregnancy Test Strip, OTC	Strip	1 test/box
FS-HCG-S30-050	Forsure hCG rapid test, strip format	Strip	50 tests/box
FS-HCG-C30-050	Forsure hCG rapid test, cassette format	Cassette	50 tests/box
FS-HCG-SH40-001b	Forsure One-Step Pregnancy Test, OTC	Strip with holder	1 test/box
HCG-S30-001	hCG rapid test, strip format *	Strip	Bulk
HCGs-S30-001	hCG rapid test (serum/urine), strip format *	Strip	Bulk
HCG-C30-001	hCG rapid test, cassette format *	Cassette	Bulk
HCG-SH40-001	hCG rapid test, strip with holder format *	Strip with holder	Bulk
HCG-M60-001	hCG rapid test, midstream format *	Midstream	Bulk
LH-C30-001	LH rapid test, cassette format *	Cassette	Bulk
LH-M60-001	LH rapid test, midstream format *	Midstream	Bulk
FSH-M60-001	FSH rapid test, midstream format *	Midstream	Bulk

* OEM/Private labels. We provide specialised private labelling consultation, registration and manufacturing services. Please speak to our sales & marketing team.



FORSURE self-test / OTC packages.



Midstream, one-step, cassette and strip formats.

TS Diagnostics Sdn Bhd

No 7 & 9, Jalan TPP 5/11,
Taman Perindustrian Puchong,
47160 Puchong, Selangor,
Malaysia.

Ph number +6 03 8068 7311
Fax number +6 03 8051 1778
contactus@forsure-test.com
www.forsure-test.com

