



# **FORSURE**

## LATEX KIT REAGENTS

FORSURE Latex Kit offer a quick and simple assay to diagnose a range of pathogens and medical conditions. The assay is based on an immunological reaction between the detected analyte in the sample and its corresponding antibody or antigen already coated on latex particles

FORSURE Latex Kit cover a selection of routine tests in serology and microbiology. They are conveniently packed in sizes of 50 or 100 tests and include all the necessary reagents, controls, stirrers and slides to conduct the test. FORSURE Latex Kit is affordable, easy to use, dependable and offer a clear and visible agglutination for doubt-free results within 2 minutes.

# FORSURE

## LATEX KIT REAGENTS

Product Code	Latex Kits	Test Method	Size
FS-ASO-L50-001	Forsure ASO Latex Kit	Slide Test	50 tests
FS-ASO-L100-001	Forsure ASO Latex Kit	Slide Test	100 tests
FS-CRP-L50-001	Forsure CRP Latex Kit	Slide Test	50 tests
FS-CRP-L100-001	Forsure CRP Latex Kit	Slide Test	100 tests
FS-DDM-L50-001	Forsure D-Dimer Latex Kit	Slide Test	50 tests
FS-DDM-L100-001	Forsure D-Dimer Latex Kit	Slide Test	100 tests
FS-IM-L50-001	Forsure IM Latex Kit	Slide Test	50 tests
FS-IM-L100-001	Forsure IM Latex Kit	Slide Test	100 tests
FS-RF-L50-001	Forsure RF Latex Kit	Slide Test	50 tests
FS-RF-L100-001	Forsure RF Latex Kit	Slide Test	100 tests
FS-SLE-L50-001	Forsure SLE Latex Kit	Slide Test	50 tests
FS-SLE-L100-001	Forsure SLE Latex Kit	Slide Test	100 tests
FS-TOXO-L50-001	Forsure Toxo Latex Kit	Slide Test	50 tests
FS-TOXO-L100-001	Forsure Toxo Latex Kit	Slide Test	100 tests
FS-STAPH-L50-001	Forsure Staphylococcus Latex Kit	Slide Test	50 tests
FS-STAPH-L100-001	Forsure Staphylococcus Latex Kit	Slide Test	100 tests
FS-STREP-L50-001	Forsure Streptococcus Latex Kit	Slide Test	50 tests
FS-STREP-L100-001	Forsure Streptococcus Latex Kit	Slide Test	100 tests



### TS Diagnostics Sdn Bhd

No 7 & 9, Jalan TPP 5/11,  
Taman Perindustrian Puchong,  
47160 Puchong, Selangor,  
Malaysia.

Ph number +6 03 8068 7311  
Fax number +6 03 8051 1778  
contactus@forsure-test.com  
www.forsure-test.com

